



BURTON

Burton SL3mz Slit Lamp



Since 1947, Burton has set the standard for exceptional quality, value, service and support.

The Burton SL3mz is a sure fire winner in design and function.

- Classic Zeiss type slit lamp
- Outstanding Galilean optics
- Moisture-proof, mildew-proof anti-reflective lens
- 0 mm to 14mm slit width
- Drum change 10X, 16X, and 25X magnifications
- Optional external illumination system
- Single-handed elevation system
- Rheostat for smooth light intensity



* Applanation Tonometer not included.

The Burton SL3mz combines traditional design with superior optics, the SL3mz will be a welcome addition to your practice.

Distributed By: