



BURTON

Burton-XL4000 Chair & Stand



Since 1947, Burton has set the standard for exceptional quality, value, service and support.

No hydraulics, means a cleaner office

Quiet, smooth electrical operation of both the chair and the slit lamp arm.

Modular slit lamp arm for ease of installation and maintenance.

Instrument Wells with Charge Indicator

Variable Overhead Lamp
Power & Intensity Control on Stand

Motorized Adjustable Instrument Arm.
Additional Arm & Chair Adjustment
Switch located on Arm (to facilitate
wheel chair patients diagnosis)

Multifunction chair & Unit Elevation Control



Movement control is powered with a quiet, German DC-Marked motor. Rotation is a full 180 degrees. Incline can be controlled from both the chair side and control panel. Installation can be for either right hand or left hand operations.

Chair is fully flexible, rotating, and reclinable. The seat is polyurethane filled for the ultimate seating pleasure. The headrest is adjustable with a single locking lever.

The foot switch has fluid movement to all positions up, down, front, back and return.



Distributed By: