



exceptional
INSTRUMENT STAND

Model 2000-ST

RIGHTmed



MODEL 2000-ST INSTRUMENT STAND

A revolutionary mix of design, function & value, the 2000-ST is the right choice for today's modern practice.

FEATURING:

modern DESIGN

all metal construction with stainless steel accents provide a modern, pleasing style.

logical INSTRUMENT CONSOLE

durable membrane switches, logically positioned, control power mains, chair up/down and auto-return, overhead lamp with rheostat, and binding post (BIO) control with rheostat. Three rechargeable wells accommodate all brands of hand held instruments.

intelligent SLIT LAMP ARM

an engineering achievement, our fully counterbalanced arm, incorporates the full function and range, with the lock/release mechanism conveniently located at your finger tips. A 10" (26cm) range of travel and counter balance range up to 60lbs (27.2kg), this arm easily accommodates all, including wheel chair patients. The chair up/down switch is conveniently located on the front of the arm and is reversible.

ergonomic UPPER (REFRACTOR) ARM

engineered with natural, and comfortable hand movements in mind, the lock/release lever engages a single lock mechanism in all three dimensions (forward/backward/vertical). A 12.5" (31.75cm) vertical range, and counter balance up to 20lbs (9.1kg) ensures precise and effortless positioning.

elegant HALOGEN OVERHEAD LIGHT

the overhead light completes the modern appearance of the 2000-ST, yet offers a long-lasting halogen lamp, providing bright, direct illumination for general lighting or focused minor procedures. Metal lamp arm construction allows stable, one handed movements with superior durability.

optional THIRD ARM

with the same ergonomically designed lever mechanism as the refractor arm, the optional third arm delivers a vertical range of 12.5" (31.75cm) and counter balance up to 30lbs (13.6kg) to easily accommodate the addition of a Keratometer to your examination room.

sensible OPTIONAL ADDITIONS



STORAGE TRAY



WHEELCHAIR KIT



BIO HANGER



Specifications: base 22.25" (56.5cm) wide x 23.5" (59.7cm) depth, height to top of console 34.5" (87.6cm) height to top of pole 68" (172.7cm)

Shipping weight: 255lbs (115.7kg) gross 255lbs (102.1kg) net. With a policy of continuous improvement, we reserve the right to make changes to specifications without notice.

Distributed by:



769 Wesley Ave.
Tarpon Springs Fl. 34689
Local: 727-934-0130
Toll Free 866-934-1030
www.Patternless.com